

MANUAL FOR THE NURSE CALL BINARY STATION

Working with the NCB																												
Flavours	<p>The NCB comes with a number of pre settable operation modes called flavours that couple the input signals to the indicators.</p> <p>Change of one or more inputs, triggers the NCB, depending on the set flavour, to set the indicators and transmit a frame with a two bytes payload.</p> <ul style="list-style-type: none"> - Byte 0 contains the code for the flavour - Byte 1 contains the status of the extender inputs <p>Also the Alarm flag in de NTM status byte (byte 11) of the data frame is set when a channel is activated and reset when no channel is activated.</p> <table border="1"> <thead> <tr> <th>Input</th> <th>Off</th> <th>Payload byte 1</th> </tr> </thead> <tbody> <tr> <td>Pull cord</td> <td>Change of state input</td> <td>0000 000X</td> </tr> <tr> <td>Emergency button</td> <td>Change of state input</td> <td>0000 00X0</td> </tr> <tr> <td>Reed switch</td> <td>Change of state input</td> <td>0000 0X00</td> </tr> <tr> <td>Extra button</td> <td>Change of state input</td> <td>0000 X000</td> </tr> <tr> <td>Aux1</td> <td>Change of state input</td> <td>000X 0000</td> </tr> <tr> <td>Aux2</td> <td>Change of state input</td> <td>00X0 0000</td> </tr> <tr> <td>Aux3</td> <td>Change of state input</td> <td>0X00 0000</td> </tr> <tr> <td>Aux4</td> <td>Change of state input</td> <td>X000 0000</td> </tr> </tbody> </table>	Input	Off	Payload byte 1	Pull cord	Change of state input	0000 000X	Emergency button	Change of state input	0000 00X0	Reed switch	Change of state input	0000 0X00	Extra button	Change of state input	0000 X000	Aux1	Change of state input	000X 0000	Aux2	Change of state input	00X0 0000	Aux3	Change of state input	0X00 0000	Aux4	Change of state input	X000 0000
Input	Off	Payload byte 1																										
Pull cord	Change of state input	0000 000X																										
Emergency button	Change of state input	0000 00X0																										
Reed switch	Change of state input	0000 0X00																										
Extra button	Change of state input	0000 X000																										
Aux1	Change of state input	000X 0000																										
Aux2	Change of state input	00X0 0000																										
Aux3	Change of state input	0X00 0000																										
Aux4	Change of state input	X000 0000																										
Offuse	<p>If the NCB is to be used as an output device and power consumption needs to be low, put the NTM into synchronized sleep mode by setting the sleep flag and actor flag. The NTM will try to find a beacon signal to synchronize with, the moment it leaves program mode (J13 removed).</p>																											
Door	<p>Flavours</p> <p>In the Offuse mode the NCB can be used as a general purpose I/O device. Channel values from expander are transmitted in byte 1. The setting of the channels is done with the AUXP command.</p> <p>Payload byte 0 = 79</p> <p>Typical application for the “Door” flavour is the door/window alarm contact. Either use a single magnet mounted directly next to the reed switch or use a separate magnet/switch set and connect it to the external switch terminals. Auxiliary channels can be set either as input or output. Inputs will trigger transmission of the state of the NCB.</p> <p>Outputs can be controlled via the CPRT command</p> <p>Payload byte 0 = 68</p>																											
Indicator																												
Alarm	<table border="1"> <thead> <tr> <th>Output</th> <th>At rest</th> <th>Activated</th> </tr> </thead> <tbody> <tr> <td>LED 1</td> <td>Off</td> <td>Off</td> </tr> <tr> <td>LED 2</td> <td>Off</td> <td>Quick flash</td> </tr> </tbody> </table>	Output	At rest	Activated	LED 1	Off	Off	LED 2	Off	Quick flash																		
Output	At rest	Activated																										
LED 1	Off	Off																										
LED 2	Off	Quick flash																										

MANUAL FOR THE NURSE CALL BINARY STATION

Working with the NCB																
<p>Nsalarm Syncalarm</p>	<p>The Indicator flavour is meant as wireless over the door side indicator. The application program must link an Indicator to a master Syncalarm NCB by providing the device ID of the master using the STMS command. The indicators will follow the indicators on the master NCB. Payload byte 0 = 73</p> <p>Alarm flavour makes the station applicable as a two button warning station. Pull cord switch and external button act parallel as button 1. Emergency button act as button 2. Button 1 signals help (twice pulled), button 2 signals emergency.</p> <p>The station is reset by activating the reed switch or input Aux1. Payload byte 0 = 65</p> <table border="1"> <thead> <tr> <th>Output</th> <th>At rest</th> <th>Button 1 once</th> <th>Button 1 twice</th> <th>Button 2</th> </tr> </thead> <tbody> <tr> <td>LED 1</td> <td>Off</td> <td>Off</td> <td>Slow flash</td> <td>Quick flash</td> </tr> <tr> <td>LED 2</td> <td>Off</td> <td>Off</td> <td>Slow flash</td> <td>Quick flash</td> </tr> </tbody> </table> <p>Flavours Nsalarm and Syncalarm perform the same function. Syncalarm NCB broadcasts a beacon signal, Nsalarm does not.</p>	Output	At rest	Button 1 once	Button 1 twice	Button 2	LED 1	Off	Off	Slow flash	Quick flash	LED 2	Off	Off	Slow flash	Quick flash
Output	At rest	Button 1 once	Button 1 twice	Button 2												
LED 1	Off	Off	Slow flash	Quick flash												
LED 2	Off	Off	Slow flash	Quick flash												
<p>Triple</p>	<p>They are similar to flavour Alarm, using a two button warning system. Pull cord switch and external button act parallel as button 1. Aux1 input acts as button 2. Alternative: - Emergency button act as emergency button</p> <p>Pushing button 1 once or more, generates a help request Pushing button 1 together with 2, generates an emergency request Pushing button 2 only, resets the alarm Alternative: - Using the emergency button resets the alarm</p> <p>Payload byte 0 = 78,83</p> <table border="1"> <thead> <tr> <th>Output</th> <th>At rest</th> <th>Button 1 once</th> <th>Byte1 twice</th> <th>Button 1&2</th> </tr> </thead> <tbody> <tr> <td>LED 1</td> <td>Off</td> <td>Slow flash</td> <td>Slow flash</td> <td>Quick Flash</td> </tr> <tr> <td>LED 2</td> <td>Off</td> <td>Off</td> <td>Slow flash</td> <td>Quick Flash</td> </tr> </tbody> </table> <p>Flavour Triple pull is a variant of the Nurse call flavours. It is about slow or fast pull cord switching. The external switch acts again parallel to the pull cord switch.</p> <p>1 pull within means help 2 pulls within 10 seconds means assistance required 3 pulls within 10 seconds is emergency 2 pulls later than 10 seconds after the last pull resets the system</p>	Output	At rest	Button 1 once	Byte1 twice	Button 1&2	LED 1	Off	Slow flash	Slow flash	Quick Flash	LED 2	Off	Off	Slow flash	Quick Flash
Output	At rest	Button 1 once	Byte1 twice	Button 1&2												
LED 1	Off	Slow flash	Slow flash	Quick Flash												
LED 2	Off	Off	Slow flash	Quick Flash												

MANUAL FOR THE NURSE CALL BINARY STATION

Working with the NCB																
	<p>Alternatives:</p> <ul style="list-style-type: none"> - Reed switch resets the alarm - Emergency button acts as emergency button <table border="1" style="width: 100%; text-align: center;"> <thead> <tr> <th>Output</th> <th>At rest</th> <th>pulled once</th> <th>Pulled twice</th> <th>Pulled trice</th> </tr> </thead> <tbody> <tr> <td>LED 1</td> <td>Off</td> <td>Slow flash</td> <td>Slow flash</td> <td>Slow Flash</td> </tr> <tr> <td>LED 2</td> <td>Off</td> <td>Off</td> <td>Slow flash</td> <td>Quick Flash</td> </tr> </tbody> </table>	Output	At rest	pulled once	Pulled twice	Pulled trice	LED 1	Off	Slow flash	Slow flash	Slow Flash	LED 2	Off	Off	Slow flash	Quick Flash
Output	At rest	pulled once	Pulled twice	Pulled trice												
LED 1	Off	Slow flash	Slow flash	Slow Flash												
LED 2	Off	Off	Slow flash	Quick Flash												
<p>Parameters/commands</p> <p>NCB programming commands</p> <p>n= number</p> <p>LF = line feed</p> <p>B= bit number 0 – 7</p> <p>X= 0 or 1</p> <p>_ = mandatory space</p>	<p>In addition the existing commands for the NTM (see application note 1), there are 9 extra commands to set parameters for the NCB.</p> <p>TEMP?LF returns actual temperature</p> <p>TMMX?LF request temperature threshold</p> <p>TMMX= nLF temperature alarm threshold</p> <p>DTdt?LF request temperature rise setting</p> <p>DTdt=nLF minimum temperature rise alarm</p> <p>NCFG?LF request NCB flavour setting</p> <p>NCFG=xLF setting of the NCB flavour</p> <p style="padding-left: 20px;">x = O: Offuse (direct extender I/O)</p> <p style="padding-left: 20px;">D: Door (door contact)</p> <p style="padding-left: 20px;">I: Indicator (side indicator)</p> <p style="padding-left: 20px;">A: Alarm (basic call buttons)</p> <p style="padding-left: 20px;">N: Nsalarm (no sync call station)</p> <p style="padding-left: 20px;">T: Tripple (triple pull call station)</p> <p style="padding-left: 20px;">E: External (external I2C device)</p> <p style="padding-left: 20px;">S: Syncalarm (beacon emitter for side indicator)</p> <p>AUXP?LF request auxiliary channel setting</p> <p>AUXP=B_XLF program the four auxiliary channels into inputs and outputs</p> <p style="padding-left: 40px;">For Aux1 choose B = 4 (not available with flavours that use Aux1)</p> <p style="padding-left: 40px;">For Aux2 choose B = 5</p> <p style="padding-left: 40px;">For Aux3 choose B = 6</p> <p style="padding-left: 40px;">For Aux4 choose B = 7</p> <p style="padding-left: 40px;">Example:</p> <p style="padding-left: 40px;">AUXP=4_1 makes Aux1 input</p> <p style="padding-left: 40px;">AUXP=4_0 makes Aux1 output</p> <p>SRMS?LF request address of master</p> <p>SRMS=nLF set address of associated master</p>															

MANUAL FOR THE NURSE CALL BINARY STATION

Working with the NCB																												
	<p>RDEX=<i>nLF</i> Read external I2C device</p> <p>WDEX=<i>nLF</i> Write external I2C device</p> <p>CPRT?<i>LF</i> reports the status of the 4 fixed inputs</p> <p>CPRT=<i>nLF</i> sets the 4 aux outputs when defined as outputs. When a channel is set to input this command is ignored</p> <p>“n” is a decimal value as indicated in following table. Sum of these values will activate multiple outputs (active low).</p> <table border="1" style="margin-left: auto; margin-right: auto;"> <thead> <tr> <th>Function</th> <th>n</th> <th>Exterder bit</th> </tr> </thead> <tbody> <tr> <td>Pull cord</td> <td>No use</td> <td style="background-color: red;">0</td> </tr> <tr> <td>Emergency button</td> <td>No use</td> <td style="background-color: red;">1</td> </tr> <tr> <td>Reed switch</td> <td>No use</td> <td style="background-color: red;">2</td> </tr> <tr> <td>Extra button</td> <td>No use</td> <td style="background-color: red;">3</td> </tr> <tr> <td>Aux1</td> <td>16</td> <td style="background-color: orange;">4</td> </tr> <tr> <td>Aux2</td> <td>32</td> <td style="background-color: green;">5</td> </tr> <tr> <td>Aux3</td> <td>64</td> <td style="background-color: green;">6</td> </tr> <tr> <td>Aux4</td> <td>128</td> <td style="background-color: green;">7</td> </tr> </tbody> </table>	Function	n	Exterder bit	Pull cord	No use	0	Emergency button	No use	1	Reed switch	No use	2	Extra button	No use	3	Aux1	16	4	Aux2	32	5	Aux3	64	6	Aux4	128	7
Function	n	Exterder bit																										
Pull cord	No use	0																										
Emergency button	No use	1																										
Reed switch	No use	2																										
Extra button	No use	3																										
Aux1	16	4																										
Aux2	32	5																										
Aux3	64	6																										
Aux4	128	7																										
Additional information	<p>Datasheet NTM_3</p> <p>Datasheet Nurse call binary station</p> <p>Application note 1; programming the NTM</p> <p>Application note 2; Ninthway high secure radio network</p>																											